

Sample Preparation for *In Situ* Enzymatic Digestion in SDS Page Gels

- Wear gloves, lab coat and work as dust free as possible.
- We recommend that gels be stained with Coomassie Blue (either R250 or G250). In the case of G250 (colloidal) Coomassie Blue, the sensitivity of staining appears to be about 25 ng (*i.e.*, 0.3 pmol) transferrin in a 0.5 mm thick gel. Currently, we do *not* recommend the use of silver staining. If you need to use silver stain contact me for a protocol.
- Scan the gel and send a copy of the image with the excised band circled and molecular weight markers labeled.
- Rinse all tubes with methanol and wipe down all surfaces with methanol. This is important to avoid keratin contamination.
- Wipe razor blades with methanol-soaked lint-free cloth
- Excise the band from the gel in such a manner as to avoid removing excess gel that does not contain any protein. 1mm cubes.
- Place the excised band in a clear Eppendorf tube and freeze. Do not dry it or leave the gel band in destain buffer or in any other liquid.
- Excise a similar size piece of "blank" gel (that does not contain any protein) and put it in a separate Eppendorf tube so that this blank section of gel can be used as a control to identify HPLC artifact peaks.
- Be sure to include as much information about the sample as you can such as species, type of preparation (nuclear extract, BALF) etc.

Updated July 2003
Charlotte Mobarak, PhD.
BMSB 249
Department of Biochemistry and Molecular Biology
University of New Mexico
505-272-8193
cmobarak@salud.unm.edu